



610.756.7879


[Home](#) » [Products](#) » [Merit Rochester-Pean Forceps](#)

MERIT ROCHESTER-PEAN FORCEPS

Product Number: 97-471

5-1/2" Merit Rochester-Pean Forceps Curved

Merit™ Rochester-Pean Hemostatic Forceps are used for clamping larger tissue and vessels for hemostasis. The full horizontal serrations and availability in multiple lengths make these forceps a versatile instrument used in multiple procedures. They are available curved or straight. Merit™ instruments are Sklar's mid-grade line of products and most commonly utilized for routine and in-office minor procedures. These cost effective, high quality instruments are crafted from Pakistan steel and come with a 10 year guarantee. Merit™ instruments deliver an economical solution for physicians and practices that require quality, durable instruments. This product is curved with full, horizontally serrated jaws, and a length of 5-1/2 inches.



Additional Information

Surname	Rochester-Pean
Length	5-1/2"
Instrument Type	Hemostatic Forceps
Curvature	Curved
Working Surface Style	Serrated
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • Service@SklarCorp.com