



610.756.7879

[Home](#) » [Products](#) » [Merit Adson-Brown Tissue Forceps, 7x7 Teeth - 4-3/4"](#)

## MERIT ADSON-BROWN TISSUE FORCEPS, 7X7 TEETH - 4-3/4"

Product Number: 97-603

4-3/4" Merit Adson-Brown Tissue Forceps, Side Grasping 7x7 Teeth

Merit™ Adson-Brown Tissue forceps are for holding and manipulating delicate tissues. These forceps have a wide, flat thumb grasp area that is commonly serrated. The jaws are short and the tips are narrow. They are similar to standard Adson forceps, however, Adson-Brown forceps feature multiple rows of teeth, typically 7x7 or 9x9. They are used to grasp delicate tissue, most commonly in plastic surgery. Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. These cost effective, high quality instruments are crafted from Pakistan steel and come with a 10 year guarantee. Merit™ instruments deliver an economical solution for physicians and practices that require quality, durable instruments. This product is straight with 7x7 teeth, and a length of 4-3/4 inches.



### Additional Information

Surname	Adson-Brown
Length	4-3/4"
Tip Configuration	7x7 Teeth Side-Grasping
Instrument Type	Tissue Forceps
Curvature	Straight
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)