



610.756.7879


[Home](#) » [Products](#) » [Merit Iris Dressing Forceps](#)

MERIT IRIS DRESSING FORCEPS

Product Number: 98-129

4" Merit Iris Dressing Forceps, Half-Curved, Serrated

Iris Dressing Forceps are thumb style forceps commonly used for placing or removing dressings in many types of procedures. Although Dressing Forceps were originally designed to specifically hold materials such as gauze, they have come to have a variety of additional functions. They may be used for removing foreign objects or dead tissue from wounds, as well as other small surgical materials such as tubing. Serrated jaws assist in having a secure hold on items. Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. These cost effective, high quality instruments are crafted from Pakistan steel and come with a 10 year guarantee. Merit™ instruments deliver an economical solution for physicians and practices that require quality, durable instruments. This product is curved with smooth, blunt tips and a length of 14 inches.



Additional Information

Surname	Iris
Length	4"
Tip Configuration	Blunt
Instrument Type	Dressing Forceps
Curvature	Slightly Curved
Working Surface Style	Serrated
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • Service@SklarCorp.com