



610.756.7879


[Home](#) » [Products](#) » [Merit Iris Tissue Forceps](#)

## MERIT IRIS TISSUE FORCEPS

Product Number: 98-146

4" Merit Iris Tissue Forceps 1x2 Teeth Half-Curved

Merit™ Iris Tissue Forceps are most commonly used for grasping and handling tissue in ophthalmic procedures. The narrow tips, however, allow this instrument to be beneficial for use in other microsurgical procedures as well. Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. These cost effective, high quality instruments are crafted from Pakistan steel and come with a 10 year guarantee. Merit™ instruments deliver an economical solution for physicians and practices that require quality, durable instruments. This product is slightly curved with 1x2 teeth and a length of 4 inches.



### Additional Information

Surname	Iris
Length	4"
Tip Configuration	1x2 Teeth
Instrument Type	Tissue Forceps
Curvature	Slightly Curved
Working Surface Style	Smooth
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)