



610.756.7879

[Home](#) » [Products](#) » [Merit Curette-Excavator](#)

MERIT CURETTE-EXCAVATOR

Product Number: 98-505

5-1/2" Merit Curette-Excavator Double End, 1mm/2mm



Additional Information

Length	5-1/2"
Tip Dimensions	1mm/2mm
Tip Configuration	Blunt
Instrument Type	Excavator Curette
Curvature	Straight
Working Surface Style	Smooth
Handle	Octagon
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • Service@SklarCorp.com