



610.756.7884

[Home](#) » [Products](#) » [Operating Scissors - 6"](#)

OPERATING SCISSORS - 6"

Product Number: 13-1060

6" Operating Scissors, Sharp/Sharp

Sklar® Operating Scissors, also known as Surgical Scissors, are a versatile surgical instrument. With three available tip configurations (sharp/sharp, blunt/blunt, or sharp/blunt), varying curvatures and multiple lengths, these scissors are used for general, multipurpose cutting and dissecting. Sharp operating scissors cut tougher tissue and are frequently used for dissecting, whereas blunt operating scissors are more often used for cutting materials such as sutures, tubing, or gauze. This product is straight with smooth, sharp/sharp tips and a length of 6 inches.



Additional Information

Length	6"
Tip Configuration	Sharp/Sharp
Instrument Type	Operating Scissors
Curvature	Straight
Working Surface Style	Smooth
Handle	Finger Rings
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com