



610.756.7884

[Home](#) » [Products](#) » [Cushing Brain Forceps](#)

## CUSHING BRAIN FORCEPS

Product Number: 22-6870

7" Cushing Brain Forceps Delicate 1x2 Teeth

Sklar® Cushing Tissue Forceps are light to medium weight thumb forceps used on delicate tissue. They are commonly used in neurosurgery. Cushing forceps have very narrow tips and are useful in small surgical areas. They feature 1 X 2 teeth to securely grasp the tissue, but also have a stop peg to prevent an overly harsh grasp that may crush the tissue. They are available with a standard, Gutch, or bayonet handle depending on the surgical site and user preference. This product is straight with delicate 1x2 teeth, and a length of 7 inches.



### Additional Information

Surname	Cushing
Length	7"
Tip Configuration	1x2 Teeth
Instrument Type	Tissue Forceps
Curvature	Straight
Working Surface Style	Blunt
Pattern	Delicate
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)